

EQUIPMENT	BRAND/MODEL	TECHNICAL/DETAILS	FACILITY - ANCON	FACILITY - BCI	FACILITY - BOCAS DEL TORO	FACILITY - COIBITA	FACILITY - FORTUNA	FACILITY - GALETA	FACILITY - GAMBOA	FACILITY - NAOS	FACILITY - TUPPER
BioSafety Cabinet, Class II	Labconco Biologic Class II	Hepa filters 99.99+%. Overall depth x height: 31.2" x 61.7" (79.2 x 156.7 cm)							2	1	1
Chemical Fume Hood	Hamilton, Concept								2	1	
Chemical Fume Hood	N/A							1			
Chemical Fume Hood	AMS (Air master systems corp.).	Fume Hood Safety, 4 ft³			2						
Deionizer	Deionizer. Bantman	Typical use of water: for Autoclaves & dishwashers. Max flow rate: 38 L/hr.							2		
Horizontal Flow clean bench	Baker EG-6252, Laminar 6"	HEPA filters 99.99+%. Provides product protection. Does not provide any personnel or environmental protection.								1	
Laminar Flow Chamber	Envirco ESM.	With HEPA filter that remove all particles 0.3 microns and larger in size with an efficiency of 99.99%		1							
Laminar Flow Chamber	SAS. Sentry Air Systems, Inc	Laminar Flow Microbiological Safety Cabinet. HEPA filters			1						
pH meter	Mettler Toledo S220	pH measuring range: -2 to 20.00. pH accuracy (±): 0.002						1	1	1	
pH meter	ThermoScientific, Orion Star A111	pH meter benchtop. pH measuring range: -2.00 to 16.00. pH accuracy (±): 0.01			1						
pH meter	ThermoScientific, Orion 3 Star	pH Portable. pH measuring range: -2.00 to 19.999. pH accuracy (±): 0.002			1						
Water Distiller	Barnstead MP 6A	Product Water capacity: 6.1 L/ hr. Four electric heaters							1		
Water Distiller	Barnstead, Fistream II	Product Water capacity: 1.4 L/hr. Two electric heaters		1			1	1			
Water Distiller	Barnstead, MP 11	Product Water capacity: 13 L/hr. Four electric heaters		1							
Water Distiller	Thermo Scientific.		1								1
Water Distiller	Barnstead.									1	
Water Purification	EMD Millipore Water purification.									1	
Water Purification	Thermo scientific, Barnstead	Ultrapure cartridge and final ultrapure filter of 0.2 µm.							2		
Water Purification	Thermo Scientific		1								
Water Purification	Barnstead Intl. Model 4642-33	E-PURE, 4 cartridge filtration. Pre-treatment, High Capacity & Ultra Pure, final 0.2 micron filter			1						

Contacts per facility					
Tupper / Fortuna / ANCON	BCI	Bocas del Toro	Galeta	Gamboa	Coibita / Naos
Raineldo Urriola	Melissa Cano	Plinio Gondola	Lourdes Vargas	Raineldo Urriola	Juan Maté
Scientific Coordinator	Scientific Coordinator	Scientific Coordinator	Scientific Coordinator Assistant	Scientific Coordinator	Manager for Scientific Affairs & Operations
Tel 212-8124 Cell 6614-9831	Tel 212-8927 Cell 6675-9248	Tel 212-8569 Cell 6611-3270	Tel 212-8194 Cell 6330-2820	Tel 212-8124 Cell 6614-9831	Tel 212-8253 Cel 6140-3847
urriolar@si.edu	canom@si.edu	Gondolap@si.edu	vargasl@si.edu	urriolar@si.edu	matej@si.edu
Felix Rodriguez				Isis Ochoa	
Scientific Coordinator Assistant				Scientific Coordinator Assistant	
Tel 2128134 Cell 67719138				Tel 2128453 Cell 67290436	
rodriguezf@si.edu				ochoai@si.edu	